

VIGILANT+labelsTM

EVERY LABEL, EVERY TIME

CLICK-TO-COMPLY SOLUTIONTM

The labeling of all anesthesia medications, medication containers, and other solutions is a risk reduction activity consistent with safe medication management. So why do the bare minimum?

The cost of implementing a labeling system is minuscule when compared to the cost associated with patient injury or death. Adverse drug events cost Americans \$2 billion each year. At least 1.5 million Americans are injured each year due to adverse drug events, of which 7,000 result in death.

WHY CHOOSE US

Vigilant Labels specializes in patient and OR safety and efficiency, and ensures that our labels are compliant with all regulatory agencies.

- Vigilant Labels Click-to-ComplyTM solution works out-of-the-box with no support needed from IT.
- Vigilant Labels Click-to-ComplyTM solution labeling matches ASA, ASTM and ISMP standards and guidelines.
- Vigilant Labels Click-to-ComplyTM solution mirrors current OR workflow, standardizing medication dosages, and saving time.

**MEDICATION ERROR IS THE MOST COMMON MEDICAL ERROR
DON'T LET THIS HAPPEN AT YOUR HOSPITAL**

OUR VALUE TO YOU

At Vigilant Labels, we work to free physicians and nurses from the mundane task of creating labels so that they can focus on what's really important - the patient.

1. Safety

Studies have shown that consistent labeling of syringes can decrease the number of medication errors caused by anesthesia providers.

2. Compliance

Anesthesia providers are more likely to use a system like the Click-to-ComplyTM solution due to its ease of use for each and every case.

3. Efficiency

The Vigilant Label system can quickly print standardized labels individualized for your facility and anesthesia providers.

Click-to-Comply™ - Physical Specifications

Dimensions:	16.8 cm x 18.0 cm x 21.5 cm (6.61 in x 7.09 in x 8.46 in)
Weight:	1.72 kg (3.80 lbs)
Power Supply Requirements:	100–240 VAC, 50–60 Hz, 60 W
Operating Temperature:	5 °C to 40 °C (41 °F to 104 °F)
Storage Temperature:	-40 °C to 60 °C (-40 °F to 140 °F)
Relative Humidity:	20 to 80% non-condensing

Click-to-Comply™ - Print Speed, Resolution, and General Media Specifications

Print Mode:	Direct thermal
Print Speed (variable):	203 dpi: 50 to 200 mm/sec (2 to 8 ips) 300 dpi: 50 to 150 mm/sec (2 to 6 ips)
Print Resolution:	203 dpi: 8 dots/mm 300 dpi: 11.8 dots/mm
Print Width:	203 dpi: 54 mm (2.1 in) max 300 dpi: 56 mm (2.2 in) max
Media Width:	19.05 to 118 mm (0.75 to 4.65 in)
Label Roll Diameter (maximum):	127 mm (5 in)
Label Length (minimum):	6.35 mm (0.25 in)
Gap Length Between Labels:	2 to 10 mm (0.08 to 0.39 in)
Black Mark Length:	3 to 10 mm (0.12 to 0.39 in) centered between tickets, 15% maximum reflectance.
Label Stop Sensor:	Adjustable from left media edge to 4.5 mm right of center on LCD printers

SERVICE PLUS PACKAGE

Vigilant Labels Service Plus Package is designed to make the solution simple to maintain. Our services package will deliver maximum uptime for anesthesia and nursing workflows virtually eliminating the need for IT Support oversight. This comprehensive service offering supports both hardware and software no matter what happens during the course of usage in the OR or on the patient floors. Selecting Service Plus delivers the ultimate in investment protection.

ABOUT US

At Vigilant Labels, we are committed to helping you improve patient safety and hospital compliance. Our elegant labeling solution can be customized to meet the needs of your patients in the clinical, perioperative, and hospital settings. Through our stand-alone label printer system, you can print a variety of different Joint Commission compliant labels for your patients.

CONTACT US

- (855) 522-3573
- <http://vigilantlabels.com/contact-us/>
- <http://vigilantlabels.com>

